

Shimadzu UK's new 'Summer Package Deal' serves up sizzling GCMS system savings

[Previous](#) | [Next](#)

This summer Shimadzu UK is offering up to 25 percent savings on a complete GCMS system that is ready to produce high quality analytical results for a diverse range of laboratory applications, which cover R&D, academia, hospitals and general chemical manufacturing.



Shimadzu UK's 'Summer GCMS Package Deal' features the Shimadzu QP2010SE GCMS plus a capillary column supported with a NIST library for a special price tag of only £30,000. Another 'Summer Deal' sees Shimadzu's QP2010 Ultra EI GCMS being offered for as little as £40,000 even though the system is fully equipped with a capillary column and benefits from NIST library support.

Both the QP2010SE and QP2010 Ultra EI are rugged GCMS systems offering reliable performance and featuring a dual filament with auto changeover which ensure minimal data is lost during sample sequencing. Extremely versatile, robust and easy to maintain, both systems enable front opening access to the source area. Simplifying maintenance still further, the systems are designed to ensure that the quadrupole analyser is not disturbed during source cleaning.

The high performance QP2010 Ultra EI, which is capable of rapidly pumping down to enable qualitative analysis within 10 to 15 minutes, also achieves industry leading fast scan rate of 20,000 amu/sec and offers the highest signal to noise performance of >500:1 for 1pg OFN.

The performance of both the QP2010SE and QP2010 Ultra EI systems can be easily enhanced with the addition of a direct sample inlet probe accessory that can heat up to 500°C and offers total software control.

Some examples of the diverse range of analytical applications that can be achieved with the QP2010SE and QP2010 Ultra EI systems include analysing semi volatiles in drinking water, identifying PBDEs in plastics and creating fatty acid profiles in neo-natal screening for inborn errors of metabolism.

To qualify for the 'Summer GCMS Package Deal' organisations need to place an order with Shimadzu UK for the Shimadzu QP2010 system or Shimadzu's QP2010 Ultra EI system before 30th September 2011.

For more information visit www.shimadzu.co.uk

Already a member? [Log in](#)

Interested? Require further information?

Note. Your details will be referred to the company and they will provide you with more information regarding your enquiry directly

If you have not logged into the website then please enter your details below.

About You

Prof Dr Mr Mrs Miss Ms

First Name

Last Name

Email

Send Information To

Organization

Organization Address

Zip / Postal Code

--Country--

Telephone Number

Job Title

Primary Specific Discipline

Work Field

Type of enquiry

Message

I am looking to purchase this product in:

One Month



Please upgrade to a [supported browser](#) to get a reCAPTCHA challenge.

[Why is this happening to me?](#)

Request Information

Related Articles:

- [Shimadzu UK grows sales orders by 30%!](#)
- [Shimadzu launches new GCMS-TQ8030 ultra fast triple quadrupole GCMS for high-end applications](#)
- [Shimadzu UK launches liquid chromatography mass spectrometer \(LCMS\) offering class-leading sensitivity and speed](#)
- [Shimadzu Training Courses: Chromatography & Mass Spectrometry](#)
- [Opportunities and challenges in the booming North American cannabis market](#)
- [5 Days of Exciting Science and Innovative Talks by Experts in the World of Liquid Phase Separation Science and Analysis](#)
- [Advances in Bio-Separations - ChromSoc release Preliminary Programme](#)
- [Industry-first glycolipid-specific library](#)
- [‘Sharing best practice’ events at Shimadzu University help participants improve their research methods](#)
- [The forefathers of Shimadzu win the 2012 Pittcon Heritage Award](#)
- [World’s fastest LC-MS/MS from Shimadzu delivers increased productivity and sensitivity performance](#)
- [Shimadzu’s compact diffraction grating innovation creates new benchmark in](#)

[molecular spectroscopy](#)

- [New corporate claim represents new mission statement](#)
- [New multimode inlet enables a single GC-MS system to support multiple analysis tasks](#)
- [New Compact, Research-grade UV-VIS Spectrophotometers Deliver Industry's Highest Absorbance Level and Ultra-Low Stray Light](#)
- [Shimadzu's New LabSolutions DB Provides Secure, Efficient Data Management for LC/GC Laboratories](#)
- [Shimadzu offers Universal Testing Machines direct to UK customers](#)
- [Shimadzu releases new version of GC/MS Method Package for Residual Pesticides in Foods](#)
- [Fully supported Edman sequencer offers cost effective N-Terminal protein sequencing benefits](#)

Newsletter Sign up

Subscribe here



Subscribe to receive our newsletters for the latest news on new laboratory products, research, Industry news and more



Weekly Update | Separation Science | Microscopy & Image Analysis | Monthly Update

Popular this Month...

Our Top 10 most popular articles this month

Today's Picks...

Looking for a Supplier?

Search by company or by product

Company Name:

Product:

S E A R C H

Please note Lab Bulletin does not sell, supply any of the products featured on this website. If you have an enquiry, please use the contact form below the article or company profile and we will send your request to the supplier so that they can contact you directly.

Lab Bulletin is published by newleaf marketing communications ltd

[Previous](#) | [Next](#)

[Back to top](#)